

## ● SAFETY TIPS

### • Matters needing attention

- (1) Although internal resistance is high due to the weight and design of the working head, continuous vibration may cause tension in the user's hand.
- (2) In order to avoid heat accumulation in the heat sink, hands or articles do not block the heat sink.
- (3) In special circumstances or accidents, please stop using it immediately and consult professionals.
- (4) Turn off the power immediately and stop using if any of the following abnormal conditions occur during use;
  - Water and foreign matter enter the body;
  - The body, power cord, accessory wire and plug have abnormal heat and leakage.
  - Abnormal voice;
  - Abnormal smell and smoke;
  - Any attachments are damaged;
  - Natural disasters such as thunder and earthquakes occur suddenly.
- (5) After use, to cut off the power supply, when not in use must be unplugged from the power socket.
- (6) It is forbidden to use the impact massager in the following environment: high temperature and humidity environment, dusty environment, oil smoke or water vapor environment.
- (7) In case of failure in use, the power should be turned off immediately, and it must be disassembled and repaired by professionals or sent to dealers and manufacturers for repair, not repaired by yourself.
- (8) The instrument is a massage instrument, not a medical instrument.

### • Taboo

- (1) Patients with heart disease, cardiac pacemaker, hypertension, acute inflammation, deep vein thrombosis, goiter, asthma, cancer, etc.;
- (2) people with blood vessel rupture and bleeding, trauma or bleeding wound;
- (3) pregnant women;
- (4) areas with plastic surgery and artificial fillers; A person with metal material in the body (including metal teeth)
- (5) skin inflammation symptoms, body edema;
- (6) Abnormal immune system function;
- (7) people with numbness and no response to heat;
- (8) Not supported in the vicinity of certain tissues, eyes and surrounding areas, myocardium, Spinal cord, gonads, kidney pulp and liver;

# Product maintenance

## 1. Check before use

Before using the impact massage machine, the following checks should be carried out:

- Check for any mechanical damage;
- Inspect all exposed wires, insert parts and accessories;
- Check the function of all possible instruments for human use and ensure that they are in good working condition.

## 2. Maintenance and maintenance of impact massager

- The impact massager is a general mechanical and electronic product with a service life of 5 years under normal working conditions (the warranty period is 1 year).
- Users should wipe and clean the dust on the surface of the massage instrument at least once a week. Wipe with a clean cloth or soft brush. Alcohol, gasoline or other detergent is prohibited. Wipe the host and physiotherapy handle, and liquid is prohibited to flow into the physiotherapy instrument.
- Frequent maintenance of working heads and accessories can improve their service life. It is recommended to maintain them once every month.
- When the massage instrument fails, it should be stopped in time and repaired by professional maintenance personnel to prevent the expansion of the fault and the occurrence of safety accidents.

Technology brings  
you beauty

# Product specification

## V101 Intelligent impact massager

Promotes cell division and improves blood circulation




## Product Overview

Electromagnetic dot array shock wave uses dc pulse current to generate a magnetic field through the coil and promote the motion of the shock film to generate shock wave, so as to realize focusing, flat emission, or scattering force deep into the fascia layer to break the adhesion of muscle and fascia and carry out corresponding conditioning.


For people with chronic muscle strain, a lattice shock wave can be painful and jarring. Because the muscle has been stuck to the fascia for a long time, the position of the muscle and fascia has become distorted, and to remove the adhesion of the muscle and fascia requires a force that can penetrate the fascia layer to remove them and return them to the correct natural position.

## Product parameters



Model: V101	Color: blue and white
Voltage: AC110V - AC220V	Frequency:50Hz/60Hz
power: 200W MAX	current: 5A MAX
Set energy level: 5-250mj	Set frequency: 1-16Hz
Strike Times, shift gears: 5	Number of jobs: 19
Net weight: 5kg	Gross weight : 6kg
Instrument size: 35*28.8*20cm	

## Working head display



## Instrument connection steps

(1) Insert the multi-core plug at the end of the handle connection wire into the OUT socket on the front of the host device.



(2) Insert the power cord plug of the host into the power socket of the host device; Turn on the main power switch at the rear of the host device.



(3) Select the right working head, hold the handle with one hand, unscrew the working head from the handle with the other hand in the (counterclockwise) direction, and then tighten the working head in the (clockwise) direction; Operate handle one button switch, press one button to open, press again to close.

## Password protection and language setting



The initial password of the instrument is **23456**. In this interface, long press the handle switch to display the password reset interface.

Operation interface can switch to four languages: simplified Chinese, traditional Chinese, English and Spanish.

### Language choice



### • ED physical therapy

- (1) The frequency of ED physiotherapy is 1-10Hz, and the recommended frequency is 3-8Hz.
- (2) Please adjust the energy according to user requirements.
- (3) Operation from low to high.

## APPLICATION SUGGESTIONS

The user can adjust the energy level and frequency according to the severity of the damage

Apply physical therapy advice					
Region	Frequency/level	Hit Count	Working hours	Total hits	Interval time
Elbow	3Hz/70mj	2,000-2,500	3-5Minutes	6,000-12,500	5-7days
Hand/wrist	1Hz/70mj	2,000	3-5Minutes	6,000-10,000	5-7days
Pubic symphysis	3Hz/90mj	2,000-2,500	5Minutes	6,000-12,500	Daily 1 day
Knee	3Hz/90mj	2,000	3-5Minutes	6,000-10,000	5-7days
Feet/ankles	3Hz/70mj	1,000	3-5Minutes	30,00-5,000	1 days
Shoulder	6Hz/120mj	2,000-2,500	3-5Minutes	6,000-12,500	5-7days
Lumbar Vertebra	8Hz/150mj	2,000-2,500	3-5Minutes	6,000-12,500	5-7days
Fat Accumulation	10Hz/150mj	2,000-2,500	5Minutes	6,000-12,500	5-7days
Buttocks	4Hz/180mj	2,000-2,500	5Minutes	6,000-12,500	5-7days